



# Surgical & Wound Care Catalogue 2023


International Sales



# Product range and areas of use

Medela products comply with international industry standards: ISO 9001 and ISO 13485 as well as the American Good Manufacturing Practice (GMP) guidelines. Our products are CE marked in accordance with the Medical Device Directive 93/42/EEC.

**Hyperlinks in the chart below guide you directly to the product section on our homepage.**



## Quick Access

### New feature

Scan or click to access the product section on our homepage offering you:

- Detailed product information and therapy insights
- Downloads/Videos
- and many more

## PVS



Surgical and Airway Suction



Vacuum Assisted Delivery



Fluid Collection

## CT



Cardiothoracic Drainage

## NPWT



Negative Pressure Wound Therapy

## Go-to Market Tool-Box



# Table of Contents

## Quick Access

Links to PVS, CT, NPWT and Market Tool-Box	4
--	---

## Cardiothoracic Drainage

Thopaz+	6
Disposables	7
Drainage catheters	10
Accessories	12

## Airway and Surgical Vacuum Suction Solution

Clario Airway Suction Pump	14
Clario Toni Airway Suction Pump	15
Accessories	16
Vario 18 Portable Suction Pump	17
Vario 8 Low Vacuum, Portable Suction Pump	18
Vario c/i Portable Suction Pumps for Constant or Intermittent Suction	19
Accessories	20
Dominant Flex Surgical Suction Pump	21
Basic Surgical Suction Pumps	22
Accessories	24
Medela Reusable Collection System – Standard Setup Recommendation	25
Medela RCS (Reusable Collection System)	26
Medela Disposable Collection System – Standard Setup Recommendation	30
Medela DCS (Disposable Collection System)	31
Accessories	32
Filters	33
Consumables	35
Patient Tubing	36
Vacuum Tubing	38

## Vacuum Assisted Delivery

Obstetric Product Range – Single Use	40
Obstetric Product Range – Reusable	40
Vacuum Extractor Pump	43
Accessories	43

## Negative Pressure Wound Therapy System (NPWT)

Invia Liberty	44
Disposables	44
Accessories	45
Invia Motion	46
Disposables	46
Accessories	46
Invia Dressing Kits	47

## Promotion Material

Promotion Material	49
--------------------	----

## Languages/Remark/Signs and Symbols

Languages	50
Remark	50
Signs and Symbols	50

# Thopaz+ Cardiothoracic Drainage

## Medela takes chest drainage therapy to a new level of care



- Compact, portable vacuum source
- Regulated pressure at the patient's chest
- Digital display of air leak, fluid and pressure
- Notifications on system status and catheter patency indicator
- Data recording and transfer to PC
- Warranty: Device 2 years
- Battery life min. 4 hours, average 10 hours
- Weight 1 kg
- WxHxD: 170x95x235 mm
- Suction performance 5L/min
- Pressure range 0 to -10 kPa



Scan or click for  
further information

**079.1000**

### **Thopaz+ – Digital cardiothoracic drainage system – Basic unit 200.5718**

**Consisting of:**

- 087.0059 Thopaz+ wall mains adapter/charger
- 079.1007 Thopaz+ CD
- 200.6841 Instructions for Use Thopaz+ EN/DE/FR/IT/NL/SV/ES/PT

**079.1002**

### **Thopaz+ – Digital cardiothoracic drainage system – Basic unit 200.5718**

**Consisting of:**

- 087.0059 Thopaz+ wall mains adapter/charger
- 079.1007 Thopaz+ CD
- 200.6842 Instructions for Use Thopaz+ EN/DA/FI/NO/CS/PL/SK/EL

**079.1003**

### **Thopaz+ – Digital cardiothoracic drainage system – Basic unit 200.5718**

**Consisting of:**

- 087.0059 Thopaz+ wall mains adapter/charger
- 079.1007 Thopaz+ CD
- 200.6843 Instructions for Use Thopaz+ EN/JA/KO/RU/TR

# Thopaz+ Disposables

## Canisters without solidifier

Compatible to Thopaz/Thopaz+

**079.0011**      **Thopaz canister 0.3L, sterile, antifoaming**      **12 units per PC**



**079.0016**      **Thopaz canister 0.8L, sterile, antifoaming**      **6 units per PC**



**079.0018**      **Thopaz canister 2.0L, sterile, antifoaming**      **12 units per PC**



## Canisters with solidifier

Compatible to Thopaz/Thopaz+

**079.0012**      **Thopaz canister 0.3 L with solidifier, sterile, antifoaming**      **12 units per PC**



**079.0017**      **Thopaz canister 0.8 L with solidifier, sterile, antifoaming**      **6 units per PC**



**079.0019**      **Thopaz canister 2.0 L with solidifier, sterile, antifoaming**      **12 units per PC**



# Thopaz+ Disposables

## Tubings with secretion extraction point







Compatible to Thopaz/Thopaz+

<b>079.0021</b> 	<b>Thopaz tubing single, sterile</b> – Compatible with 20–32Fr* – Length 150 cm	<b>10 units per PC</b> 
<b>079.0022</b> 	<b>Thopaz tubing double, sterile</b> – Compatible with 20–32Fr* – Length 150 cm	<b>10 units per PC</b> 
<b>079.0023</b> 	<b>Thopaz tubing single, small connector, sterile</b> – Compatible with 12–20 Fr* – Length 150 cm	<b>10 units per PC</b> 
<b>079.0024</b> 	<b>Thopaz tubing double, small connector, sterile</b> – Compatible with 12–20 Fr* – Length 150 cm	<b>10 units per PC</b> 
<b>079.0025</b> 	<b>Thopaz tubing single, large connector, sterile</b> – Compatible with 28–36 Fr* – Length 150 cm	<b>10 units per PC</b> 
<b>079.0026</b> 	<b>Thopaz tubing double, large connector, sterile</b> – Compatible with 28–36 Fr* – Length 150 cm	<b>10 units per PC</b> 

\*Information depends on the catheter design. If you have any questions, please contact your local sales representative.



# Thopaz + Disposables

101039326	<b>Tubing extension, 1 connector, sterile</b>	<b>15 units per PC</b>
	<b>The set includes:</b> <ul style="list-style-type: none"><li>- 1 silicone tube (45 cm)</li><li>- 1 straight connector</li></ul>	
101039327	<b>Tubing extension, 3 connectors</b>	<b>10 units per PC</b>
	<b>The set includes:</b> <ul style="list-style-type: none"><li>- 3 silicone tubes (30 cm)</li><li>- 3 straight connectors</li><li>- 1 Y-connector</li></ul>	
079.0039	<b>Thopaz sealing cap, sterile</b>	<b>20 units per PC</b>
	<ul style="list-style-type: none"><li>- For closing unused connectors</li><li>- Compatible with all tubing</li></ul>	

## Drainage catheters

100008ATS (8 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100010ATS (10 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100012ATS (12 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100014ATS (14 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100016ATS (16 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100020ATS (20 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100024ATS (24 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100028ATS (28 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100032ATS (32 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC
100032ATS (32 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , straight	10 units per PC



- Made of silicone
- Minimises the risk of clot formation due to the innovative ClotStop<sup>®</sup> coating made of hydrogel polymer (PVP)
- Scaled
- Sizes: 8–12 Fr, length 22.9 cm  
14–36 Fr, length 50.8 cm

111016ATS (16 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , right angle	10 units per PC
111120ATS (20 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , right angle	10 units per PC
111124ATS (24 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , right angle	10 units per PC
111128ATS (28 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , right angle	10 units per PC
111132ATS (32 Fr)	Thoracic catheter with ClotStop <sup>®</sup> , right angle	10 units per PC



- Made of silicone
- Minimises the risk of clot formation due to the innovative ClotStop<sup>®</sup> coating made of hydrogel polymer (PVP)
- Scaled
- Sizes: 16–32 Fr, length 50.8 cm

111075AT	Mediastinal drain with ClotStop <sup>®</sup> , (7 mm, eyes on sides)	10 units per PC
111095AT	Mediastinal drain with ClotStop <sup>®</sup> , (9 mm, eyes on sides)	10 units per PC
111099AT	Mediastinal drain with ClotStop <sup>®</sup> , (9 mm, eyes on top and bottom)	10 units per PC





- Made of silicone
- Specifically developed for cardiac surgery
- Flattened oval shape for maximum drainage flow with the minimum catheter profile
- Minimises the risk of clot formation due to the innovative ClotStop<sup>®</sup> coating made of hydrogel polymer (PVP)
- Sizes: 7 mm (O.D. 12.57 x 6.81 mm; I.D. 8.00 x 5.81 mm), length 50.8 cm and  
9 mm (O.D. 16.51 x 9.41 mm; I.D. 13.36 x 8.64 mm), length 50.8 cm
- Versions: Eyes on sides or eyes on top and bottom

Exclusively distributed in Europe by Medela

 Axiom Medical Inc.  
19320 Van Ness Ave.  
Torrance, CA 90501  
www.axiommed.com



# Drainage catheters

111016S (16 Fr)	<b>Thoracic catheter, straight</b>	10 units per PC
100020S (20 Fr)	<b>Thoracic catheter, straight</b>	10 units per PC
100024S (24 Fr)	<b>Thoracic catheter, straight</b>	10 units per PC
100028S (28 Fr)	<b>Thoracic catheter, straight</b>	10 units per PC
	<ul style="list-style-type: none"> <li>- Made of silicone</li> <li>- Scaled</li> <li>- Sizes: 16–28 Fr, length 50.8 cm</li> </ul>	
111016S (16 Fr)	<b>Thoracic catheter, curved</b>	10 units per PC
111120S (20 Fr)	<b>Thoracic catheter, curved</b>	10 units per PC
111124S (24 Fr)	<b>Thoracic catheter, curved</b>	10 units per PC
111128S (28 Fr)	<b>Thoracic catheter, curved</b>	10 units per PC
	<ul style="list-style-type: none"> <li>- Made of silicone</li> <li>- Scaled</li> <li>- Sizes: 16–28 Fr, length 50.8 cm</li> </ul>	

Exclusively distributed in Europe by Medela

 Axiom Medical Inc.  
19320 Van Ness Ave.  
Torrance, CA 90501  
[www.axiommed.com](http://www.axiommed.com)

 0413

# Thopaz+ Accessories

087.0059

## Thopaz+ wall mains adapter/charger



– With different plugs (Euro-GB-US-AUS/NZ)

*Compatible to Thopaz/Thopaz+*

079.1006

## Thopaz+ carrying strap



– For carrying Thopaz/Thopaz+

*Compatible to Thopaz/Thopaz+*

079.0037

## Docking station – Basic unit 200.2233



– For storing Thopaz+ and charging the battery  
– Can be attached to standard rails, beds\* and IV stands\*  
\*when using a holder with a standard rail

*Compatible to Thopaz/Thopaz+*

079.0038

## Adapter to docking station for Thopaz canister 2L



– Recommended for attaching the 2L canister

*Compatible to Thopaz/Thopaz+*

079.0034

## USB cable



– For connecting Thopaz+ to the PC  
–Type: 2.0 A/M- Mini B/M

*Compatible to Thopaz/Thopaz+*

079.0036

## Holder with standard rail



079.0040

## Universal holder with standard rail



– Length standard rail 30cm

# Thopaz + Accessories

077.1456

**Vacuum gauge – 10 kPa for Thopaz safety test**

– Pressure range 0 to –10 kPa



## Training material

079.0033

**Siphon Demo Tool**

# Clario Airway Suction Pump



- 3 Vacuum settings
  - min = -135 mmHg/-18 kPa
  - med = -270 mmHg/-36 kPa
  - max = -600 mmHg/-80 kPa
- Suction capacity 15l/min.
- SafetyChamber for absolute safety
- Canister set included
- Threefold hygienic design
- Membrane system ensures no overflow into pump
- Warranty: Device 2 years / battery 6 months



Scan or click for further information

## Clario, Standard version

- AC pump
- Instructions for use

**Disposable finger tubing Art. no. 074.0003 (15 units per PC) must be ordered separately (not included within the pump set).**

Numbers	Instructions for use	Voltage	Plug
014.0010	177.1405 (2x)*: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0011	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0014	177.1406: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240 V/50-60 Hz	Euro
014.0040	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	UK
014.0050	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	AUS/NZ
014.0020	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	US

\* only for Austria/Germany/Switzerland

## Clario, AC/DC version

- AC/DC pump
- Instructions for use

**Disposable finger tubing Art. no. 074.0003 (15 units per PC) must be ordered separately (not included within the pump set).**

Numbers	Instructions for use	Voltage	Plug
014.0110	177.1405 (2x)*: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0111	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0114	177.1406: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240 V/50-60 Hz	Euro
014.0120	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	US
014.0140	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	UK
014.0150	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	AUS/NZ

\* only for Austria/Germany/Switzerland

## Clario, AC/DC version, Deluxe

- AC/DC pump
- 1 spare canister set
- SoftCase carrying bag
- Instructions for use

**Disposable finger tubing Art. no. 074.0003 (15 units per PC) must be ordered separately (not included within the pump set).**

Numbers	Instructions for use	Voltage	Plug
014.0210	177.1405 (2x)*: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0211	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	Euro
014.0214	177.1406: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240 V/50-60 Hz	Euro
014.0240 (incl. car cable)	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	UK
014.0250	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	AUS/NZ
014.0230	177.1405: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240 V/50-60 Hz	US

\* only for Austria/Germany/Switzerland

# Clario Toni Airway Suction Pump



- Vacuum settings
  - min = -70 mmHg/-9 kPa for newborns
  - med = -115 mmHg/-15 kPa for children
  - max = -600 mmHg/-80 kPa for viscous secretions
- Suction capacity 15 l/min.
- SafetyChamber for absolute safety
- Canister set included
- Threefold hygienic design
- Membrane system ensures no overflow into pump
- Warranty: Device 2 years/battery 6 months



Scan or click for further information

## Clario Toni, AC/DC version

- AC/DC pump
- Instructions for use

Disposable finger tubing Art. no. 074.0003 (15 units per PC) must be ordered separately (not included within the pump set).

Numbers	Instructions for use	Voltage	Plug
014.0022	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	AUS/NZ
014.0023	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	Euro
014.0024	200.8401: PL/DA/EL/CS/SL/HU/ET/LV/LT/SL	100-240 V/50-60 Hz	Euro
014.0025	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	UK
014.0032	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	US

## Clario Toni, AC/DC version, Deluxe

- AC/DC pump
- 1 spare canister set
- SoftCase carrying bag
- Instructions for use

Disposable finger tubing Art. no. 074.0003 (15 units per PC) must be ordered separately (not included within the pump set).

Numbers	Instructions for use	Voltage	Plug
014.0026	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	AUS/NZ
014.0027	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	Euro
014.0028	200.8401: PL/DA/EL/CS/SL/HU/ET/LV/LT/SL	100-240 V/50-60 Hz	Euro
014.0029 (incl. car cable)	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	UK
014.0031	200.8400: EN/DE/FR/IT/NL/SV/NO/FI/PT/ES	100-240 V/50-60 Hz	US

# Clario and Clario Toni Accessories


074.0001

## Canister set



### Consisting of:

- Suction canister, 0.5 litre
- Lid with plugs
- Floater

 (Minimum order quantity: 10 units)




074.0002

## SafetyChamber



### Consisting of:

- SafetyChamber housing
- Lid to SafetyChamber
- 2 valve discs

 (Minimum order quantity: 10 units)




014.0030

## SafetyChamber Toni



### Consisting of:

- SafetyChamber housing
- Lid to SafetyChamber
- 2 valve discs

 (Minimum order quantity: 10 units)



074.0003

## Disposable fingertip tubing, sterile

15 units per PC



### Consisting of:

- Length 150 cm
- CH 25



074.0017

## Car connection cable

- For operating the Clario AC/DC pump in the car



074.0007

## SoftCase carrying bag

- For the safe transport of the home care system with accessories
- With 2 carrying straps for attaching to wheelchairs or as a backpack and carrying handle





# Vario 18 Portable Suction Pump



- Max. vacuum up to -75 kPa (-563 mmHg)
- Suction capacity 18l/min.
- QuatroFlex™ technology
- Membrane vacuum regulator, vacuum gauge in kPa and mmHg
- Integrated holder on one side for additional suction jar
- Including 5 disposable filters
- Silicone tubing Ø 7 x 12 mm, 25 cm with 1 coupling piece and 1 coupling piece 90°
- Instruction for use
- Warranty: Device 2 years/QuatroFlex 6 months
- WxHxD: 170x380x285 mm

Medela recommends a 1 litre reusable or 1.5 litre disposable collection system



Scan or click for further information

## Vario 18 AC

Numbers	Instructions for use	Voltage	Plug
026.2110	177.2604 (2x)*: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	230-240V/50Hz	Euro
026.2110	177.2605 (1x): ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	230-240V/50Hz	Euro
026.2112	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	230-240V/50Hz	Euro
026.2112	177.2605: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	230-240V/50Hz	Euro
026.2150	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	230-240V/50Hz	AUS/NZ
026.2120	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	120 V/60Hz	US
026.2160	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT 101040029: KO	230-240V/60Hz	Euro

\* only for Germany/Austria/Switzerland



- Same as above but battery version
- Can be operated from mains or battery mode (approx. 30 min.) or with 12V car-connection
- Battery is charged during mains operation
- Optical and acoustic warning signal when battery low
- Warranty: Device 2 years/Battery 6 months/QuatroFlex 6 months
- WxHxD: 170x380x285 mm

Medela recommends a 1 litre reusable or 1.5 litre disposable collection system

## Vario 18 AC/DC

Numbers	Instructions for use	Voltage	Plug
026.5110	177.2604 (2x)*: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60Hz	Euro
026.5110	177.2605 (1x): ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240V/50-60Hz	Euro
026.5112	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT 101040029: KO	100-240V/50-60Hz	Euro
026.5112	177.2605: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240V/50-60Hz	Euro
026.5140	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60Hz	UK
026.5120	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60Hz	US

\* only for Germany/Austria/Switzerland

# Vario 8 Low Vacuum, Portable Suction Pump



- Vacuum range from -10 to -100 mbar (-1 to -10 kPa)
- Suction capacity 8l/min.
- QuatroFlex™ technology
- Membrane vacuum regulator, vacuum gauge in kPa and mbar
- Integrated holder on one side for additional suction jar
- Silicone tubing Ø 7 x 12 mm, 25 cm with 1 coupling piece and 1 coupling piece 90°
- Including 5 disposable overflow protection/bacteria filters
- Instruction for use
- Warranty: Device 2 years/QuatroFlex 6 months

Medela recommends a 1 litre reusable collection system



Scan or click for further information

## Vario 8 AC

Numbers	Instructions for use	Voltage	Plug
026.7160	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT 101040029: KO	230 – 240 V/60 Hz	Euro

\* only for Germany/Austria/Switzerland



- Same as above, but battery version
- Can be operated from mains or battery mode (approx. 120 min.) or with 12V car-connection
- Optical and acoustic warning signal when battery low
- Battery is charged during mains operation
- Warranty: Device 2 years/Battery 6 months/QuatroFlex 6 months

Medela recommends a 1 litre reusable collection system

## Vario 8 AC/DC

Numbers	Instructions for use	Voltage	Plug
026.8110	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100–240V/50–60 Hz	Euro
026.8110	177.2605: PL/DA/EL/CS/SL/HU/ET/LV/LT/SK 101040029: KO	100–240V/50–60 Hz	Euro
026.8150	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100–240V/50–60 Hz	AUS/NZ

\* only for Germany/Austria/Switzerland

# Vario c/i Portable Suction Pumps for constant or intermittent suction



- Max. Vacuum up to -55 kPa/-413 mmHg
- Suction capacity 18l/min.
- Constant and intermittent suction
- Intermittent settings ON/OFF (sec.): 8/4; 16/8; 32/16
- QuatroFlex™ technology
- Membrane vacuum regulator, vacuum gauge in kPa and mmHg
- Can be operated from mains or battery mode or with 12V car-connection
- Optical and acoustic warning signal when battery low
- Battery is charged during mains operation
- Integrated holder on one side for additional suction jar
- Including 5 overflow protection/bacteria filters
- Silicone tubing Ø 7 x 12 mm, 30 cm with 2 coupling pieces (1 with 90° angle)
- Instruction for use
- Warranty: Device 2 years/Battery 6 months/QuatroFlex 6 months

Medela recommends a 1 litre reusable or 1.5 litre disposable collection system

## Vario 18 c/i AC/DC medium vacuum

Numbers	Instructions for use	Voltage	Plug
026.8152	200.0436: EN/DE/FR/IT/SV/NO/ES/PT	100-240V/50-60Hz	AUS/NZ
026.8122	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60Hz	US



- Vacuum range from -10 to -100mbar (-1 to -10kPa)
- Suction capacity 8l/min.
- Constant and intermittent suction
- Intermittent settings ON/OFF (sec.): 8/4; 16/8; 32/16
- QuatroFlex™ technology
- Membrane vacuum regulator, vacuum gauge in kPa and mmHg
- Can be operated from mains or battery mode or with 12V car-connection
- Optical and acoustic warning signal when battery low
- Battery is charged during mains operation
- Integrated holder on one side for additional suction jar
- Including 5 overflow protection/bacteria filters
- Silicone tubing Ø 7 x 12 mm, 25 cm with 2 coupling pieces (1 with 90° angle)
- Instruction for use
- Warranty: Device 2 years/Battery 6 months/QuatroFlex 6 months

Medela recommends a 1 litre reusable collection system


## Vario 8 c/i AC/DC low vacuum

Numbers	Instructions for use	Voltage	Plug
026.8111	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT 101040029: KO	100-240V/50-60 Hz	Euro
026.8111	177.2605: ET/HU/CS/LT/LV/PL/SL/SK/DA/EL	100-240V/50-60 Hz	Euro
026.8141	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60 Hz	UK
026.8151	177.2604: EN/DE/FR/IT/NL/SV/FI/NO/ES/PT	100-240V/50-60 Hz	AUS/NZ

# Vario Accessories

077.0834

## Vario disposable filter set with overflow protection/bacteria filter



- Disposable suction tubing (50 cm) with coupling piece
- Overflow protection/bacteria filter



(Minimum order quantity: 50 units)

077.0715

## Reusable 1 L set with PSU suction jar



### Consisting of:

- 077.0110 Suction jar, polysulfone, 1 litre
- 077.0420 Large lid with conical patient connection Ø 6–10 mm and overflow protection device



**Sterile tubing Art. no. 077.0951 is not part of the set and needs to be ordered separately.**

077.0082

## Suction jar PC, 1.5 L



(Recommended order quantity: 30 units)



077.0083

## Disposable suction liner 1.5L

40 units per PC



### Packaging unit:

- Per pallet 1'440 pcs
- Per container (20 ft) 21'600 pcs
- Per container (40 ft) 43'200 pcs



(Minimum order quantity: 120 units; Standard pallet qty: 36 PCs)

077.0084

## Disposable suction liner with solidifier 1.5L

30 units per PC



### Packaging unit:

- Per pallet 1'080 pcs
- Per container (20 ft) 16'200 pcs
- Per container (40 ft) 32'400 pcs



(Minimum order quantity: 90 pieces per unit; Standard pallet qty: 36 PCs)

077.0105

## Adapter for Vario

5 units per PC



Required when using the disposable collection system with Vario



077.0831

## Vario carrying bag



077.0832

## Car connection cable for 12V DC



# Dominant Flex Surgical Suction Pump



**Dominant Flex mobile**





**Dominant Flex rack**



**Dominant Flex portable**

- Max. vacuum -95 kPa (-713 mmHg)
- Suction capacity adaptable: 40, 50, 60l/min.
- Piston/cylinder technology
- Patented flat spring
- CleanTouch on/off and flow change buttons
- Membrane vacuum regulator
- Vacuum gauge in kPa and mmHg
- Equipotential conductor on back of pump
- Safety device with silicone tubing Ø 7 x 12 mm, 60 cm with 2 coupling pieces
- Port for optional foot on/off switch (077.0723)
- Instruction for use
- Warranty: 5 years

 Portable version can be equipped with Medela suction jars 1 or 2 litres (reusable collection system) and 1.5 or 2.5 litres (disposable collection system)

 Mobile version can be equipped with all Medela suction jar sizes, reusable and disposable system



Scan or click for further information

**071.0003**

**Dominant Flex rack**

Dimensions: HxWxD: 21 x30x38 cm

**071.0002**

**Dominant Flex portable**

With carrying handle for up to 2 secretion jars

**071.0003 + 071.0034**

**Dominant Flex mobile (Dominant Flex rack + trolley)**

Numbers	Instructions for use	Countries	Plug
<b>200.5790</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT, JA	Canda, Colombia, Costa Rica, Dom. Republic, Ecuador, Japan, Mexico, Panama, Peru, Philippines, Saudi Arabia, Taiwan, Thailand, Venezuela	B
<b>200.6323</b>	US	USA	B
<b>200.5791</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Egypt, India, Jordan, Nepal, Bangladesh	D
<b>200.5651</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Austria, Belarus, Belgium, Croatia, Egypt, Finland, France, Germany, Indonesia, Iran, Italy, Jordan, Lebanon, Luxembourg, Macedonia, Montenegro, Netherlands, Norway, Pakistan, Portugal, Romania, Russia, Serbia, Spain, Syria, Sweden, Tunisia, Turkey, Vietnam	E + F
<b>200.5793</b>	ET, HU, CZ, LT, LV, PL, SK, SL DK, EL	Bulgaria, Cyprus, Czech Republic, Estonia, Greece, Hungary, Latvia, Lithuania, Poland, Slovakia, Slovenia	E + F
<b>200.5795</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Bahrain, Hong Kong, Iraq, Jordan, Kuwait, Malaysia, Oman, Malta, Pakistan, Qatar, Saudi Arabia, Singapore, Thailand, UAE, UK, Yemen	G
<b>200.5796</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Israel, Palestine	H
<b>200.5797</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Argentina, Australia, New Zealand	I
<b>200.5798</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Jordan, Switzerland	J
<b>200.5799</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Chile, Italy, Syria, Uruguay	L
<b>200.5800</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	India, Nepal, South Africa	M
<b>200.5801</b>	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Brazil	N
<b>200.5802</b>	ET, HU, CZ, LT, LV, PL, SK, SL DK, EL	Denmark	K
<b>101040188</b>	KO	Korea	E + F

# Basic Surgical Suction Pump



Basic mobile





Basic rack



Basic portable

- Max. vacuum -90 kPa (-675 mmHg)
- Suction capacity adaptable: 30l/min.
- Piston/cylinder technology
- Patented flat spring
- CleanTouch on/off buttons
- Membrane vacuum regulator
- Vacuum gauge in kPa and mmHg
- Equipotential conductor on back of pump
- Safety device with silicone tubing Ø 7 x 12 mm, 60 cm with 2 coupling pieces
- Port for optional foot on/off switch (077.0723)
- Instruction for use
- Warranty: 5 years

 Portable version can be equipped with Medela suction jars 1 or 2 litres (reusable collection system) and 1.5 or 2.5 litres (disposable collection system)

 Mobile version can be equipped with all Medela suction jar sizes, reusable and disposable system



Scan or click for further information

071.0001

**Basic rack**

Dimensions: HxWxD: 21 x 30 x 38 cm

071.0000

**Basic portable**

With carrying handle for up to 2 secretion jars

071.0001 + 071.0034

**Basic mobile (Basic rack + trolley)**

Numbers	Instructions for use	Countries	Plug
200.5803	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT, JA	Canada, Colombia, Costa Rica, Dom. Republic, Ecuador, Japan, Mexico, Panama, Peru, Philippines, Saudi Arabia, Taiwan, Thailand, Venezuela	B
200.6321	US	USA	B
200.5804	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Egypt, India, Jordan, Nepal, Bangladesh	D
200.5805	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Austria, Belarus, Belgium, Croatia, Egypt, Finland, France, Germany, Indonesia, Iran, Italy, Jordan, Lebanon, Luxembourg, Macedonia, Montenegro, Netherlands, Norway, Pakistan, Portugal, Romania, Russia, Serbia, Spain, Syria, Sweden, Tunisia, Turkey, Vietnam	E + F
200.5806	ET, HU, CZ, LT, LV, PL, SK, SL DK, EL	Bulgaria, Cyprus, Czech Republic, Estonia, Greece, Hungary, Latvia, Lithuania, Poland, Slovakia, Slovenia	E + F
200.5809	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Bahrain, Hong Kong, Iraq, Jordan, Kuwait, Malaysia, Oman, Malta, Pakistan, Qatar, Saudi Arabia, Singapore, Thailand, UAE, UK, Yemen	G
200.5810	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Israel, Palestine	H
200.5811	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Argentina, Australia, New Zealand	I
200.5812	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Jordan, Switzerland	J
200.5813	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Chile, Italy, Syria, Uruguay	L
200.5814	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	India, Nepal, South Africa	M
200.5815	DE, EN, FR, IT, NL, SV, FI, NO, ES, PT	Brazil	N
200.5816	ET, HU, CZ, LT, LV, PL, SK, SL DK, EL	Denmark	K
101040187	KO	Korea	E + F


# Dominant Flex and Basic Accessories

071.0034

## Trolley for Basic and Dominant Flex



- With standard rail 10x25 mm
- Integrated foot on/off switch
- With 4 antistatic castors, with brakes
- Cable holder included
- With equipotential conductor
- Dimensions: 51 x 41 x 77 cm
- Warranty: 1 year

 Note: not compatible with Basic 30/Dominant 50

077.0723

## Foot on/off switch



- Cord length 350 cm

077.0731

## Foot vacuum regulator



- With hook and clampholder, plastic 077.0521

077.0711

## Safety device plugged into pump



### Consisting of:

- 077.0125 Suction jar, polysulfone, 0.25 litre
- 077.0450 Lid for safety jar with overflow protection device
- 077.0922 Silicone Tubing Ø 7 x 12 mm, 60 cm with 2 coupling pieces



077.0651

## Safety device for standard rails



### Consisting of:

- 077.0125 Suction jar, polysulfone, 0.25 litre
- 077.0440 Small lid with conical patient connection, Ø 6–10 mm and overflow protection device
- 077.0521 Clampholder, plastic
- 077.0912 Silicone vacuum tubing, Ø 7 x 12 mm, 60 cm, with 1 coupling piece
- 077.0913 Silicone vacuum tubing Ø 7 x 12 mm, 100 cm with 1 coupling piece



077.0751

## Tubing holder left side



077.0752

## Tubing holder right side



# Dominant Flex and Basic Accessories

077.0188

## Trolley, language combination A



- For 5 suction jars
- Dimensions: 80x60x60 cm
- IFU: EN DE FR IT NL SV NO FI PT ES

077.0521

## Clamholder, plastic



 (Recommended order quantity: 5 pieces)

077.0523

## Wall bracket, single



 (Recommended order quantity: 20 pieces)

077.0192

## Wall adapter



077.0561

## Change-over valve

077.0563

## Change-over valve with clamholder and tubing set



### Consisting of:

- 077.0561 Change-over valve
- 077.0521 Clamholder, plastic
- 077.0921 2x Silicone Tubing Ø 7 x 12 mm, 30 cm with 2 coupling pieces

077.0562

## Disposable specimen cup for reusable system Ø 6–10 mm

10 units per PC



077.0564

## Disposable specimen cup for reusable system Ø 10–14 mm

10 units per PC





# Medela Reusable Collection System – Standard Setup Recommendation

## Setup 1 Rack, single jar



## Setup 2 Rail with change-over valve



## Setup 3 Trolley with serial connection



Scan or click for further information

	Minimal ordering Qty	Setup 1 Qty required	Setup 2 Qty required	Setup 3 Qty required
--	----------------------	----------------------	----------------------	----------------------

### Reusable Collection System (Medela RCS)

077.0130	Suction jar, polysulfone, 3L Alternatives: Suction jar, polysulfone, 0.25L, 0.5L, 1L, 2L, 5L	1 unit	1 unit	2 units	5 units
077.0420	Large lid with conical patient connection 6–10 mm, overflow protection device Alternatives: Large lid 10–14 mm or small lid 6–10 mm	1 unit	1 unit	2 units	5 units

### Patient tubing/Filter/Specimen Cup

077.0951	Disposable PVC tubing, 180 cm, Ø 12x 18 mm for patient connection Ø 5.5x8.5 mm Alternatives: Patient tubing with/without fingertip or different connectors/diameters	15 pcs units	1 unit per use	1 unit per use	1 unit per use
077.0573	Disposable bacteria filter Alternatives: Bacteria, particle or odor filters with Medela or conical connectors	units	1 unit per use	1 unit per use	1 unit per use
077.0562	Disposable specimen cup for reusable System 6–10 mm	10 units	1 unit	1 unit	1 unit

### Mount

077.0521	Clampholder, plastic	1 unit		2 units	
077.0188	Trolley	1 unit			1 unit

### Large Volume Setups

077.0563	Change-over valve with clampholder and tubing set	1 unit		1 unit	
077.0542	Reusable serial connector to connect maximum four reusable 5L jars				4 units

# Medela RCS (Reusable Collection System)

## Reusable Sets

077.0701

### Reusable 1 L set with PSU suction jar



**Consisting of:**

- 077.0110 Suction jar, polysulfone, 1 litre
- 077.0420 Large lid with conical patient connection  
Ø 6–10 mm and overflow protection device
- 077.0521 Clampholder, plastic



077.0702

### Reusable 2 L set with PSU suction jar



**Consisting of:**

- 077.0120 Suction jar, polysulfone, 2 litre
- 077.0420 Large lid with conical patient connection Ø 6–10 mm and overflow protection device
- 077.0521 Clampholder, plastic



077.0703

### Reusable 3 L set with PSU suction jar



**Consisting of:**

- 077.0130 Suction jar, polysulfone, 3 litre
- 077.0420 Large lid with conical patient connection  
Ø 6–10 mm and overflow protection device
- 077.0521 Clampholder, plastic



077.0705

### Reusable 5 L set with PSU suction jar



**Consisting of:**

- 077.0150 Suction jar, polysulfone, 5 litre
- 077.0430 Large lid with conical patient connection Ø 10–14 mm and overflow protection device
- 077.0521 Clampholder, plastic



**For Art. no. 077.0701– 077.0707, sterile tubing Art. no. 077.0951 is not part of the set and needs to be ordered separately.**

# Medela RCS (Reusable Collection System)

## Double Reusable Sets

077.0707

### Double-reusable 2 L set with PSU suction jar



**Consisting of:**

- 077.0120 2x Suction jar, polysulfone, 2 litre
- 077.0420 2x Large lid with conical patient connection  
Ø 6–10 mm and overflow protection device
- 077.0563 Change-over valve with clampholder and tubing set
- 077.0521 2x Clampholder, plastic



077.0704

### Double-reusable 3 L set with PSU suction jar



**Consisting of:**

- 077.0130 2x Suction jar, polysulfone, 3 litre
- 077.0420 2x Large lid with conical patient connection  
Ø 6–10 mm and overflow protection device
- 077.0563 Change-over valve with clampholder and tubing set
- 077.0521 2x Clampholder, plastic



077.0706

### Double-reusable 5 L set with PSU suction jar



**Consisting of:**

- 077.0150 2x Suction jar, polysulfone, 5 litre
- 077.0430 2x Large lid with conical patient connection  
Ø 10–14 mm and overflow protection device
- 077.0563 Change-over valve with clampholder and tubing set
- 077.0521 2x Clampholder, plastic




**For Art. no. 077.0701– 077.0707, sterile tubing Art. no. 077.0951 is not part of the set and needs to be ordered separately.**

# Medela RCS (Reusable Collection System)

077.0125

**Suction jar, polysulfone, 0.25 litre**



 (Recommended order quantity: 63 units)



077.0155

**Suction jar, polysulfone, 0.5 litre**




 (Recommended order quantity: 27 units)



077.0440

**Small lid with conical patient connection  
Ø 6–10 mm and overflow protection device**




 (Recommended order quantity: 36 units)



077.0450

**Lid for safety jar with overflow protection device**




 (Recommended order quantity: 36 units)



077.0110

**Suction jar, polysulfone, 1 L**




 (Recommended order quantity: 16 units)



077.0120

**Suction jar, polysulfone, 2 L**




 (Recommended order quantity: 10 units)



# Medela RCS (Reusable Collection System)

077.0130

**Suction jar, polysulfone, 3 litre**

 (Recommended order quantity: 12 units)



077.0150

**Suction jar, polysulfone, 5 litre**

 (Recommended order quantity: 9 units)



077.0420


**Large lid with conical patient connection  
Ø 6–10 mm and overflow protection device**

 (Recommended order quantity: 32 units)



077.0430

**Large lid with conical patient connection  
Ø 10–14 mm and overflow protection device**

 (Recommended order quantity: 32 units)



# Medela Disposable Collection System – Standard Setup Recommendation

## Setup 1 Rack, single jar



## Setup 2 Rail with change-over valve



## Setup 3 Trolley with serial connection



Scan or click for further information

Minimal ordering Qty

Setup 1 Qty required

Setup 2 Qty required

Setup 3 Qty required

### Disposable Collection System (Medela DCS)

077.0085	Suction jar PC, 2.5L Alternative: 077.0082 Suction jar PC, 1.5L	20 units	1	2	5
077.0087	Disposable suction liner 2.5L with solidifier Alternatives	90 units	1 unit per use	2 units per use	5 units per use
	077.0083 1.5L liner without solidifier				
	077.0084 1.5L liner with solidifier				
	077.0086 2.5L liner without solidifier				

### Patient tubing/Specimen Cup

077.0951	Disposable PVC tubing, 180 cm, Ø 12.7 x 18 mm for patient connection Ø 5.5 x 8.5 mm Alternatives: Patient tubing with/ without finger-tip or different connectors/diameters	15 units	1 unit per use	1 unit per use	1 unit per use
077.0562	Disposable specimen cup for reusable System 6–10 mm	10 units	1 unit	1 unit	1 unit

### Mount

077.0521	Clampholder, plastic	1 unit		2 units	
077.0188	Trolley	1 unit			1 unit

### Large Volume Setups

077.0563	Change-over valve with clampholder and tubing set	1 unit		1 unit	
077.0158	Serial connection patient tubing	100 units			4 units per use
077.0095	Serial connection vacuum tubing	1 unit			4 units

# Medela DCS (Disposable Collection System)

077.0082

**Suction jar PC, 1.5 L**



(Recommended order quantity: 30 units)



077.0083

**Disposable suction liner 1.5 L**

**40 units per PC**



**Packaging unit:**

- Per pallet 1'440 pcs
- Per container (20 ft) 21'600 pcs
- Per container (40 ft) 43'200 pcs

(Minimum order quantity: 120 units)



077.0084

**Disposable suction liner with solidifier 1.5 L**

**30 units per PC**



**Packaging unit:**

- Per pallet 1'080 pcs
- Per container (20 ft) 16'200 pcs
- Per container (40 ft) 32'400 pcs

(Minimum order quantity: 90 units)



077.0085

**Suction jar PC, 2.5 L**



(Recommended order quantity: 20 units)



077.0086

**Disposable suction liner 2.5 L**

**40 units per PC**



**Packaging unit:**

- Per pallet 1'440 pcs
- Per container (20 ft) 21'600 pcs
- Per container (40 ft) 43'200 pcs

(Minimum order quantity: 120 units)



077.0087

**Disposable suction liner with solidifier 2.5 L**

**30 units per PC**


















**Packaging unit:**

- Per pallet 1'080 pcs
- Per container (20 ft) 16'200 pcs
- Per container (40 ft) 32'400 pcs

(Minimum order quantity: 90 units)




# Accessories

077.0094	Specimen cup for disposable collection system	10 units per PC
		
077.0097	Y-piece Ø 7 mm/Ø 5 mm	10 units per PC
		
077.0098	Y-piece Ø 7 mm	10 units per PC
		
077.0105	Adapter for Vario	5 units per PC
	<p> Required when using the disposable collection system with Vario</p>	
077.0100	Angle piece yellow	10 units per PC
		
077.0101	Angle piece orange	10 units per PC
		
077.0102	Locking clasp	10 units per PC
		












# Filters



Use filters to protect the pump from overflow, particles and the environment from bacteria, viruses and odours. See drawing for setup.



Scan or click for further information













<p><b>101035263</b></p> 	<p><b>Disposable virus filter with overflow protection, Medela connections</b></p> <p> (Recommended order quantity: 20 PCs)</p>	<p><b>10 units per PC</b></p> 
<p><b>101035264</b></p> 	<p><b>Disposable virus filter with overflow protection, conical connections</b></p> <p> (Recommended order quantity: 20 PCs)</p>	<p><b>10 units per PC</b></p> 
<p><b>101035265</b></p> 	<p><b>Disposable virus filter, Medela connections</b></p> <p> (Recommended order quantity: 20 PCs)</p>	<p><b>10 units per PC</b></p> 

# Filters

## Bacteria Filters

<b>077.0571</b> 	<b>Disposable bacteria filter with overflow protection, Medela connections</b>  (Recommended order quantity: 20 PCs)	<b>10 units per PC</b> 
<b>077.0572</b> 	<b>Disposable bacteria filter with overflow protection, conical connections</b>  (Recommended order quantity: 20 PCs)	<b>10 units per PC</b> 
<b>077.0573</b> 	<b>Disposable bacteria filter, Medela connections</b>  (Recommended order quantity: 20 PCs)	<b>10 units per PC</b> 
<b>077.0575</b> 	<b>Disposable bacteria/odour filter, Medela connections</b>  (Recommended order quantity: 20 PCs)	<b>10 units per PC</b> 

# Consumables

077.1018	Coupling piece	10 units per PC
		
077.1019	Coupling piece 90°	10 units per PC
		
077.1017	Double conical coupling piece	10 units per PC
		
077.1020	Y-Piece	10 units per PC
		
077.1013	Clip for lid	10 units per PC
		
077.1021	Overflow protection device	10 units per PC
		

# Patient Tubing

## Disposable PVC Tubing, Ø 6.9 x 10 mm, Sterile

077.0170

150 cm

15 units per PC



- For use with Bird and Silc Cups
- Without coupling pieces

## Disposable PVC Tubing, Ø 5.5 x 8.5 mm

077.0951

180 cm, sterile

15 units per PC



- For use with Vario, Basic and Dominant Flex pumps
- Stepped conical connector and fingertip

077.0184

150 cm, sterile

15 units per PC



- For use with Vario, Basic and Dominant Flex pumps
- 2 universal funnels

077.0952

180 cm, non-sterile



- For use with Basic and Dominant Flex pumps
- 1 stepped conical connector

## Silicone Tubing, Ø 6.5 x 11.5 mm

077.0053

32 cm



- For Malmström Suction Cups
- Without coupling pieces
- For use with Vario, Basic and Dominant Flex pumps

077.0055

120 cm



- For use with Vario, Basic, and Dominant Flex pumps
- Without coupling pieces

# Patient Tubing

## Silicone Tubing, Ø 7 x 12 mm

077.0056

150 cm



- For use with Vario, Basic, and Dominant Flex pumps
- Without coupling pieces

077.0052

By the meter



- For use with Vario, Basic, and Dominant Flex pumps
- Without coupling pieces

077.0961

180 cm, green



- For use with Vario, Basic and Dominant Flex pumps
- 1 coupling pieces

077.0054

50 cm



- For use with Vario, Basic and Dominant Flex pumps
- Without coupling pieces

## Tubing for Serial Connection, Non-Sterile

077.0542

Silicone Ø 6.5 x 11.5 mm, reusable serial connector



- For use with the Reusable Fluid Collection System
- Serial connector reusable

077.0158

Disposable PVC serial connection, Ø 8 x 12 mm, patient



- For use with the Disposable Fluid Collection System 100 pcs
- Serial connector disposable

077.0160

PVC serial connection, Ø 8 x 12 mm, patient



- For use with the Disposable Fluid Collection System 1 pc
- Serial connector disposable

## Additional Tubing

077.0970










Silicone, Ø 5 x 10 mm, 25 m, non-sterile



- For use with Basic and Dominant Flex pumps
- Without coupling pieces

# Vacuum Tubing

## Silicone Tubing, Ø 7 x 12 mm, Autoclavable

<b>077.0154</b>	<b>25 cm</b>		<ul style="list-style-type: none"><li>- For use with Vario, Basic and Dominant Flex pumps</li><li>- 1 coupling piece</li></ul>
<b>077.0911</b>	<b>33 cm</b>		<ul style="list-style-type: none"><li>- For use with Vario, Basic and Dominant Flex pumps</li><li>- 1 coupling piece</li></ul>
<b>077.0912</b>	<b>60 cm</b>		<ul style="list-style-type: none"><li>- For use with Basic and Dominant Flex pumps</li><li>- 1 coupling piece</li></ul>
<b>077.0913</b>	<b>100 cm</b>		<ul style="list-style-type: none"><li>- For use with Basic and Dominant Flex pumps</li><li>- 1 coupling piece</li></ul>
<b>077.0931</b>	<b>25 cm</b>		<ul style="list-style-type: none"><li>- For use with Vario, Basic and Dominant Flex pumps</li><li>- 1 coupling piece and 1 at 90°</li></ul>
<b>077.0185</b>	<b>60 cm</b>		<ul style="list-style-type: none"><li>- For use with Vario, Basic and Dominant Flex pumps</li><li>- 1 coupling piece and 1 at 90°</li></ul>
<b>077.0921</b>	<b>30 cm</b>		<ul style="list-style-type: none"><li>- For use with Vario, Basic and Dominant Flex pumps</li><li>- 2 coupling pieces</li></ul>
<b>077.0922</b>	<b>60 cm</b>		<ul style="list-style-type: none"><li>- For use with Basic and Dominant Flex pumps</li><li>- 2 coupling pieces</li></ul>
<b>077.0189</b>	<b>300 cm</b>		<ul style="list-style-type: none"><li>- For use with Basic and Dominant Flex pumps</li><li>- 2 coupling pieces</li></ul>

# Vacuum Tubing

## Silicone Tubing, Ø 7 x 12 mm, Autoclavable

**077.0901**

**100 cm**



- For use with Basic and Dominant Flex pumps
- Without coupling pieces

**077.0902**

**200 cm**



- For use with Basic and Dominant Flex pumps
- Without coupling pieces

**077.0905**

**500 cm**



- For use with Basic and Dominant Flex pumps
- Without coupling pieces

**077.0900**

**25 m**



- For use with Basic and Dominant Flex pumps
- Without coupling pieces

**077.0054**

**50 cm**



- For use with manual vacuum extractor pump
- Without coupling pieces

**077.0095**

**25 cm**



- For serial connection of the Disposable Fluid Collection System
- With y-piece

# Vacuum Assisted Delivery / Obstetric Product Range – Single-Use



Scan or click for further information

## Obstetric Product Range – Single-Use

<p><b>077.0791</b></p>	<p><b>Single-Use Silc Cup, Ø 60 mm, sterile (EO)</b></p> <p><b>Including:</b></p> <ul style="list-style-type: none"> <li>– Plastic connector</li> <li>– With PVC tubing Ø 5.5x8.5 mm, 1.5 m</li> </ul> <p> (Order quantity: 25 Units)</p>	<p><b>5 units per PC</b></p>
<p><b>077.0792</b></p>	<p><b>Single-Use Bird Cup, Ø 50 mm, sterile (EO)</b></p> <p><b>Including:</b></p> <ul style="list-style-type: none"> <li>– Plastic connector</li> <li>– With PVC tubing Ø 5.0x7.0 mm, 1.5 m</li> </ul> <p> (Order quantity: 25 Units)</p>	<p><b>10 units per PC</b></p>

# Vacuum Assisted Delivery / Obstetric Product Range – Reusable

## Obstetric Product Range – Reusable

<p><b>077.0761</b></p>	<p><b>Silc Cup, Ø 50 mm, blue</b></p> <p><b>Including:</b></p> <ul style="list-style-type: none"> <li>– Silicone tubing Ø 6.5x11.5 mm, 1.5 m</li> <li>– 1 coupling piece</li> </ul>	
<p><b>077.0762</b></p>	<p><b>Silc Cup, Ø 60 mm, yellow</b></p> <p><b>Including:</b></p> <ul style="list-style-type: none"> <li>– Silicone tubing Ø 6.5x11.5 mm, 1.5 m</li> <li>– 1 coupling piece</li> </ul>	
<p><b>077.0078</b></p>	<p><b>Silc Cup with trumpet valve, Ø 60 mm</b></p> <p><b>Including:</b></p> <ul style="list-style-type: none"> <li>– Silicone tubing Ø 6.5x11.5 mm, 1.5 m</li> <li>– Stainless steel connector</li> </ul>	



# Vacuum Assisted Delivery / Obstetric Product Range – Reusable

## Malmström Suction Cups

077.0176

### Malmström Cup Ø 50 mm, complete



**Including:**

- Bottom-plate of stainless steel with fixed chain
- Traction handle
- Silicone tubing Ø 6.5x 11.5 mm, 32 cm and 120 cm



077.0175

### Malmström Cup Ø 60 mm, complete



**Including:**

- Bottom-plate of stainless steel with fixed chain
- Traction handle
- Silicone tubing Ø 6.5x 11.5 mm, 32 cm and 120 cm



## Components

077.0029

### Malmström Suction Cup, Ø 40 mm



**Including:**

- Bottom-plate of stainless steel with fixed chain
- Traction handle
- Silicone tubing Ø 6.5x 11.5 mm, 32 cm



077.0030

### Malmström Suction Cup, Ø 50 mm



**Including:**

- Bottom-plate of stainless steel with fixed chain
- Traction handle
- Silicone tubing Ø 6.5x 11.5 mm, 32 cm



077.0031

### Malmström Suction Cup, Ø 60 mm



**Including:**

- Bottom-plate of stainless steel with fixed chain
- Traction handle
- Silicone tubing Ø 6.5x 11.5 mm, 32 cm



# Vacuum Assisted Delivery / Obstetric Product Range – Reusable

## Bird Suction Cups

077.0173

### Bird Anterior Cup Ø 40 mm, complete



**Including:**

- Traction handle and silicone tubing Ø 6.5x 11.5 mm, 150 cm
- Bottom-plates and locking ring



077.0172

### Bird Anterior Cup Ø 50 mm, complete



**Including:**

- Traction handle and silicone tubing Ø 6.5x 11.5 mm, 150 cm
- Bottom-plates and locking ring



077.0171

### Bird Anterior Cup Ø 60 mm, complete



**Including:**

- Traction handle and silicone tubing Ø 6.5x 11.5 mm, 150 cm
- Bottom-plates and locking ring



077.0174

### Bird Posterior Cup Ø 50 mm, complete



**Including:**

- Traction handle and silicone tubing Ø 6.5x 11.5 mm, 150 cm



077.1044

### Traction handle with hook to Bird Suction Cups



## CaesAid

077.0080

### CaesAid Cup, Ø 60 mm



**Including:**

- Silicone tubing Ø 6.5x 11.5 mm, 2 m
- Stainless steel connector



# Vacuum Extractor Pump

077.0081

## Manual Vacuum Extractor Pump



- 0-90 kPa
- Consisting of:**
- One plastic vacuum bottle, 500 ml
- Vacuum gauge
- Plastic coated basket
- Silicone tubing, Ø 6.5x 11.5 mm, 50 cm and 150 cm
- Warranty: 2 years

077.0026

## Manual Vacuum Extractor Pump with Malmström Cups



- 0-90 kPa
- Consisting of:**
- Manual Vacuum Extractor Pump (see Ref. 077.0081)
- 3 Malmström Suction Cups 40, 50, 60 mm, incl. Bottom-plates of stainless steel and traction handle
- Spare gasket and membrane



077.0058

## Manual Vacuum Extractor Pump with Bird Cups



- 0-90 kPa
- Consisting of:**
- Manual Vacuum Extractor Pump (see Ref. 077.0081)
- 3 Bird Anterior Suction Cups, 40, 50, 60 mm
- 5 Bottom-plates to each anterior cup
- 5 locking rings for tubing
- 1 Traction handle
- Spare gasket and membrane



077.0060

## Manual Vacuum Extractor Pump with Bird and Silc Cups



- 0-90 kPa
- Consisting of:**
- Manual Vacuum Extractor Pump (see Ref. 077.0081)
- 3 Bird Anterior Suction Cups 40, 50, 60 mm
- 2 Silc Cups, 50 and 60 mm
- 2 Traction handles with hook
- 5 Locking rings for tubing
- 10 bottom-plates to each anterior cup
- 3m extra silicone tubing
- Spare gasket and membrane



## Accessories

017.0015

## Appar Timer for standard rail (10x25 mm) battery operated



- Liquid crystal display showing minutes and seconds (18 mm)
- Acoustic signal after 1, 5 and 10 minutes
- Recording of birth time
- Two batteries 1,5V, change of batteries without interference
- Solid grey housing (L 150 mm, W 85 mm, H 45 mm)
- 200.2336 Instructions for use, EN/DE/FR/IT/NL/SV/DK/FI/NO/ES/PT/EL
- Warranty: 2 years (except for the batteries)



# Invia Liberty Negative Pressure Wound Therapy

Available in select markets only

087.0003



- Intelligent pressure control and dynamic exudate removal from the pump to the wound site
- Adjustable pressure options from -40 to -200 mmHg, including standard -125 mmHg
- On-screen air leak indicator
- Custom air leak modes for normal and large wounds
- Easy and secure tubing attachment with the Quick-connector
- One-touch canister release button
- Choice of 0.3L and 0.8L disposable canisters with solidifier
- Battery life: average 14 hours
- Modes: constant and intermittent
- 1 mains adapter/charger
- 2 instruction manuals (1 for clinician, 1 for patient)
- Weight: 1 kg
- WxHxD: 170x95x150mm
- Warranty: 2 years



Scan or click for further information

## Invia Liberty Disposables

087.0012

**Canister with solidifier, 0.3L, sterile**

10 units per PC



- Hydrophilic overflow protection/bacteria filter with activated carbon for odor protection
- Safety chamber for filter protection against premature wetting

087.0017

**Canister with solidifier, 0.8L, sterile**

6 units per PC



- Hydrophilic overflow protection/bacteria filter with activated carbon for odor protection
- Safety chamber for filter protection against premature wetting

087.0025

**Tubing with Quick-connector**

10 units per PC



- Quick, easy and secure tubing attachment with the Quick-connector
- Length of tubing: approx. 150 cm

# Invia Liberty Accessories

079.0037

## Docking station

1 unit per PC

- W x H x D: 210 x 85 x 130 mm



079.0036

## Holder with standard rail

1 unit per PC

- Length: approx. 20.5 cm



087.0031

## Carrying bag

1 unit per PC

- WxHxD: 215x160x95 mm



087.0059

## Liberty/Thopaz+/Thopaz mains adapter/charger

1 unit per PC

- For charging the battery  
- With different plugs (Euro-GB-US-AUS/NZ)



101034390

## Plug attachment for Liberty/Thopaz+ mains adapter/charger

1 unit per PC

- Euro plug attachment



# Invia Motion Negative Pressure Wound Therapy

Available in select markets only

**087.4024 (60 days) / 087.4022 (15 days)**



- Max. vacuum: -75 mmHg
- Min. vacuum: -40 mmHg
- Battery life: average 10 hours
- Modes: constant and intermittent
- 1 mains adapter/charger
- 2 instruction manuals (1 for clinician, 1 for patient)
- 1 carrying case
- Weight: 350g
- WxHxD: 90x90x50mm
- Warranty: 2 years



Scan or click for further information

## Invia Motion Disposables

**087.4016**

**Canister/Tubing set 0.15L**

**5 units per PC**



- With solidifier, sterile
- Fast canister change with one movement
- Usable canister capacity 0.15L
- Length of tubing: approx. 150 cm

## Invia Motion Accessories

**087.4003**

**Carrying case**

**1 unit per PC**



- WxHxD: 120x210x60mm

**101033647**

**Motion mains adapter/charger**








**1 unit per PC**



- For charging the battery

# Invia Dressing Kits

Available in select markets only

<b>087.6221</b>	<b>Foam dressing kit with FitPad, Size S</b>	<b>3 units per PC</b>
<b>087.6225</b>	<b>Foam dressing kit with FitPad, Size S</b>	<b>15 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 foam pad (WxHxD: 10x8x3 cm)</li> <li>  Pore size: 600–800 microns (µm)</li> <li>  Tensile strength: 138 kPa</li> <li>- 1 transparent film (WxH: 32x26 cm)</li> <li>  Water vapour permeability: 385 g/m<sup>2</sup>/day</li> <li>- 1 FitPad</li> </ul>	
<b>087.6222</b>	<b>Foam dressing kit with FitPad, Size M</b>	<b>3 units per PC</b>
<b>087.6226</b>	<b>Foam dressing kit with FitPad, Size M</b>	<b>15 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 foam pad (WxHxD: 19x12.5x3 cm)</li> <li>  Pore size: 600–800 microns (µm)</li> <li>  Tensile strength: 138 kPa</li> <li>- 2 transparent films (WxH: 32x26 cm)</li> <li>  Water vapour permeability: 385 g/m<sup>2</sup>/day</li> <li>- 1 FitPad</li> </ul>	
<b>087.6223</b>	<b>Foam dressing kit with FitPad, Size L</b>	<b>3 units per PC</b>
<b>087.6227</b>	<b>Foam dressing kit with FitPad, Size L</b>	<b>15 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 foam pad (WxHxD: 25x15x3 cm)</li> <li>  Pore size: 600–800 microns (µm)</li> <li>  Tensile strength: 138 kPa</li> <li>- 2 transparent films (WxH: 32x26 cm)</li> <li>  Water vapour permeability: 385 g/m<sup>2</sup>/day</li> <li>- 1 FitPad</li> </ul>	
<b>087.6224</b>	<b>Foam dressing kit with FitPad, Size XL</b>	<b>3 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 foam pad (WxHxD: 60x30x1.5 cm)</li> <li>  Pore size: 600–800 microns (µm)</li> <li>  Tensile strength: 138 kPa</li> <li>- 5 transparent films (WxH: 32x26 cm)</li> <li>  Water vapour permeability: 385 g/m<sup>2</sup>/day</li> <li>- 1 FitPad</li> </ul>	
<b>087.6250</b>	<b>Abdominal dressing kit with FitPad</b>	<b>5 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 organ contact layer (oval size 60x80 cm)</li> <li>- 1 foam pad (oval size 38x25x2 cm):</li> <li>  Pore size: 600–800 microns (µm)</li> <li>  Tensile strength: 138 kPa</li> <li>- 4 transparent films (WxH: 32x26 cm)</li> <li>  Water vapour permeability: 385 g/m<sup>2</sup>/day</li> <li>- 1 FitPad</li> </ul>	
<b>087.7170</b>	<b>Gauze dressing kits with FitPad, Size M</b>	<b>3 units per PC</b>
<b>087.7172</b>	<b>Gauze dressing kits with FitPad, Size M</b>	<b>15 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 Kerlix AMD antimicrobial gauze (WxH: 17 x 16 cm)</li> <li>- 1 transparent film (WxH: 32 x 26 cm)</li> <li>- 1 FitPad</li> <li>- 1 sterile saline (30 ml)</li> </ul>	
<b>087.7171</b>	<b>Gauze dressing kits with FitPad, Size L</b>	<b>3 units per PC</b>
	<ul style="list-style-type: none"> <li>- 1 Kerlix AMD antimicrobial gauze (WxH: 370 x 11.4 cm)</li> <li>- 2 transparent films (WxH: 32 x 26 cm)</li> <li>- 1 FitPad</li> <li>- 2 sterile salines (2x 30 ml)</li> </ul>	

# Invia Dressing Kits

**087.6303**      **White foam, size S**      **10 units per PC**

– PVA foam dressing (10x7.5x1 cm)




**087.6304**      **White foam, size L**      **10 units per PC**

– PVA foam dressing (15x10x1 cm)




**087.7001**      **Invia-Silverlon**      **10 units per PC**

– Antimicrobial wound contact dressing (10x12 cm)



**087.7002**      **Invia-Silverlon**      **10 units per PC**

– Antimicrobial wound contact dressing (12x20 cm)



**087.0030**      **Invia Transparent Film**      **20 units per PC**

– WxH: 32x26 cm



**087.0028**      **Invia FitPad, sterile**      **10 units per PC**

– Length of tubing: approx. 80 cm



**087.0023**      **Invia Y-connector with Quick-connector**      **10 units per PC**

– For connecting 2 FitPads to one device





# Promotion Material

## General Promotion Material Medela

### 190.3016 Medela Pen with blue ink

(Minimum order quantity: 25 units)



### 190.3054 Medela Pocket Knife



### 190.3091 Mini Chocolates



**Packaging unit:**

– (1 carton = 370 pcs.;5.18 kg)

### 200.6959 Medela USB Stick (8 GB)



### 200.6529 Dominant Flex USB stick (4 GB)



### 200.7075 Medela Lanyard Healthcare



### 200.2552 Medela Pocket Umbrella



### 200.3017 Keyholder Cowbell



### 200.3682 Medela Water Bottle




# Languages


<b>AR</b>	Arabic	<b>EL</b>	Greek	<b>PT</b>	Portuguese
<b>BG</b>	Bulgarian	<b>HE</b>	Hebrew	<b>RO</b>	Romanian
<b>ZH</b>	Chinese simplified	<b>HU</b>	Hungarian	<b>RU</b>	Russian
<b>HR</b>	Croatian	<b>IS</b>	Icelandic	<b>SR</b>	Serbian
<b>CS</b>	Czech	<b>ID</b>	Indonesian	<b>SK</b>	Slovak
<b>DA</b>	Danish	<b>IT</b>	Italian	<b>SL</b>	Slovenian
<b>NL</b>	Dutch	<b>JA</b>	Japanese	<b>ES</b>	Spanish
<b>EN</b>	English	<b>KO</b>	Korean	<b>SV</b>	Swedish
<b>ET</b>	Estonian	<b>LV</b>	Latvian	<b>TW</b>	Taiwanese (Chinese traditional)
<b>FI</b>	Finnish	<b>LT</b>	Lithuanian	<b>TR</b>	Turkish
<b>FR</b>	French	<b>NO</b>	Norwegian		
<b>DE</b>	German	<b>PL</b>	Polish		


## Remark

Not all products are available for all markets. For verification of the availability please contact the respective sales organization or the legal manufacturer.

## Signs and Symbols

 Caution: This is a single-use product not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

 Autoclavable, for example at 134°C/10 minutes

 Not autoclavable.  
Cleaning in a cleaning and disinfecting machine

 Remarks

# Notes

A series of horizontal dotted lines for taking notes.



Medela AG  
Lättichstrasse 4b  
6340 Baar, Switzerland  
[www.medela.com](http://www.medela.com)

**International Sales**

Switzerland  
Lättichstrasse 4b  
6340 Baar  
Switzerland  
Phone +41 (0)41 562 51 51  
[customerservice.international@medela.com](mailto:customerservice.international@medela.com)  
[www.medela.com](http://www.medela.com)